



HTH Engineering, Inc.
Tarpon Springs, FL 34688
Tel: 866-254-0714 Fax: (810) 821-6940
Web: <http://www.startstop.com>

QUICK START GUIDE for:



Start-Stop[®] Pedal Power for Dragon[®] Medical One


<https://www.startstop.com/start-stop-pedal-power-for-dragon-medical-one/>

- Basic steps:**
- 1) Plug the pedal into an available USB port
 - 2) Open Dragon Medical One (or DMPE/360) and Configure HotKeys
 - 3) Enjoy and use through any connection, network or version of Dragon!
-

1) Plug the pedal into any available USB port

Once plugged in please wait approximately one minute initially for the software drivers to register from Windows services - this is an automatic process.

2) Open your Dragon software - and go to the hot keys section. For Dragon Medical One, it's under OPTIONS... Hot Keys and will look like this:



For the circled objects, set them up as shown for CONTROL (CTRL) + F11, CONTROL (CTRL) + F10 and CONTROL (CTRL) +F12

This means the pedal functions will be: Left Pedal = **Previous Field** Center Pedal = **Microphone On/Off** (On when pressed, Off when released), and Right Pedal = **Next Field**

Any hot key controllable functions of Dragon (Dragon Medical One, Dragon Medical Practice, Dragon 360) can be controlled with the pedal. The “Previous Field”, “Mic On/Off” and “Next Field” labels are silk-screened on the pedal, but any Dragon hotkey controllable functions can be implemented.

Here is the Pedal Keystrokes that are sent when the pedals are pressed:



Technical support please contact us at www.startstop.com
toll-free 866-254-0714 option 1, techsupport@startstop.com
