



# Nuance PowerMic III

**Button Configuration** – most buttons are programmable and functions can be changed.

Non-Programmable Buttons:

1. Microphone
2. Mouse Pad
3. Mouse Click Left
4. Mouse Click Right
5. Speaker
6. Enter (On back not shown)

Programmable Buttons:

1. Dictation Box
2. Previous Tab
3. Next Tab
4. Dictate
5. Previous Field
6. Accept Defaults
7. Next Field
8. Left Programmable
9. Transfer Text
10. Right Programmable

(8,9,10 Color/shape shows location.)

